

## AUTOCLAVE STERILIZER LOG

MONTH \_\_\_\_\_ YEAR \_\_\_\_\_

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
<b>DAILY</b>																																
• Clean External Surface																																
• Clean Door Gasket with Damp Cloth																																
<b>WEEKLY</b>																																
• Clean Chamber & Trays																																
• Drain Water from Reservoir & Refill																																
<b>MONTHLY</b>																																
• Clean Chamber & Plumbing																																
• Remove & Clean Filters																																
• Perform Pressure Relief Valve Check																																

INITIALS	NAME